



Recombinant Rat Serum paraoxonase/arylesterase 1 (Pon1)

Product Code	CSB-BP345837RA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55159
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	AKLLGLTLV GLVLALYKNH RSSYQTRLNA FREVTPVDLP NCTLVKGIEA GAEDLEILPN GLTFFSTGLK YPGIKSFDPS KPGKILLMDL NEKEPAVSEL AIMGNTLDMS SFNPHGISTF IDEENTVYLL VVSHPDSSST VEVFKFQEEE RSLHLKTIT HELLP SINDI AAVGPESFYA TNDHYFADPY LRSWEMYLGL SWSNVVYYSP DKVRVVADGF DFANGIGISL DGKYVYIAEL LAHKIHVYEK HANWTLTPLK VLSFDTLVDN ISVDPVTGDL WVGCHPNGMR IFFYDSENPP GSEVLRIQSI LSEDPKVTVV YAENGTVLQG TTVAAYKKGK LLIGTVFHRA LCCDL
Source	Baculovirus
Target Names	Pon1
Protein Names	Recommended name: Serum paraoxonase/arylesterase 1 Short name= PON 1 EC= 3.1.1.2 EC= 3.1.1.81 EC= 3.1.8.1 Alternative name(s): Aromatic esterase 1 Short name= A-esterase 1 Serum aryldialkylphosphatase 1
Expression Region	2-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The enzyme encoded by this gene is an arylesterase that mainly hydrolyzes paroxon to produce p-nitrophenol. Paroxon is an organophosphorus anticholinesterase compound that is produced in vivo by oxidation of the insecticide parathion. Polymorphisms in this gene are a risk factor in coronary artery disease. The gene is found in a cluster of three related paraoxonase genes at 7q21.3.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.