



Recombinant *Xenopus laevis* Transcription factor Sox-3-A (sox3-a)

Product Code	CSB-EP345910XBE-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55863
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MYSMLDTDIK SPVQQSNAPI GGPATPGGKG NASTLDQDRV KRPMNAFMVW SRGQRRKMAQ ENPKMHNSEI SKRLGADWKL LSDSDKRPF DEAKRLRAVH MKDYPDYKYR PRRKTKLLK KDKYSLPGLN LAPGVSPVAS SVGVGQRIDT YAHMNGWTNG AYSLMQDQLG YSQHPAMNSP QMQQIQHRYD MSGLQYNPMM TSAQNAYMNA AASTYSMSPA YNQQSSTVMS LASMGSVVKSEPPPPAIT SHTQRACLGD LRDMISMYP PGGDASDPSL QNSRLHSVHQ HYQSAAGPGV NGTVPLTHI
Source	E.coli
Target Names	sox3-a
Protein Names	Recommended name: Transcription factor Sox-3-A Short name= xSox3 Alternative name(s): xSox-B1
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.