



Recombinant Aspartic protease 3 (asp-3)

Product Code	CSB-MP345922CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55956
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	AFNEG LSDYSNAQYY GPVTIGTPPQ NFQVLFDGTGS SNLWVPCANC PFGDIACRMH NRFDCCKSSS CTATGASFEI QYGTGSMKGT VDNDVVCFGH DTTYCTDKNQ GLACATSEPG ITFVAAKFDG IFGMGWDTIS VNKISQPMDQ IFANSAICKN QLFAFWLSRD ANDITNGGEI TLCDDTPNHV VGNIAWEPLV SEDIWRIKLA SVVIDGTTYT SGPIDSIVDT GTSLLTGPTD VIKKIQHKIG GIPLFNGEYE VECSKIPSLP NITFNLGGQN FDLQGKDYIL QMSNGNGGST CLSGFMGMDI PAPAGPLWIL GDVFIGRFYS VFDHGNKRVG FATSRTGK
Source	Mammalian cell
Target Names	asp-3
Protein Names	Recommended name: Aspartic protease 3 EC= 3.4.23.-
Expression Region	56-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.