



Recombinant Escherichia coli Uncharacterized protein yfjL (yfjL)

Product Code	CSB-YP345971ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52127
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MESNDSGGVA AKHGFLFQDC VAAYHVTRML RDKTIRSVRC EVTDDIDIVS DGYIDFVQVK STGKTRWNIS DIVQNSKGAD KKTIPCSSL HKSMQCESDL SLGRRYSIVT EEKVNKTLEY LTISPNARLD KPGRQELIDD LNKRTDNFLT DSGISVSDWI DAATWEVFSS LRELELLGIK NIRLASQDLH GVILSSETVA EDIWCRILDV VTRKGEHSRR IHSADDKSYL RPDLLEWFKQ RVEDDQSRSG RKIYVKRDLP HILTPFRAPM ASVCAKRKGQ VLHQQYSLKK YRYKHIADNV CQWLDEVFLR PKEMSDIHLK TFIKRRERLK NSVFKSLHDV SEFLGRVLLH ATIRQHHESQ PIPCMLYVEK AGAEKILENV HIVRRDPEGD QLWIGFSELV TDINIAVRLP EIRDQLYEDI SDCIDTARKK ILDIKDDNYL LRHDIDEILD GSQPFDAHLD RFTFVLFVGY DSNLLTEPET PGFEDDLEKE TAVLFEKFAA DLIEDSPFAN LCIHVFIYPA PSLERLTQLV DEKVREVV
Source	Yeast
Target Names	yfjL
Protein Names	Recommended name: Uncharacterized protein yfjL
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.