



Recombinant Arabidopsis thaliana 3-oxoacyl-[acyl-carrier-protein] synthase I, chloroplastic (KAS1)

Product Code	CSB-BP345988DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52410
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	ASTV SAPKRETPDK KRVVITGMGL VSVFGNDVDA YYEKLLSGES GISLIDRFDA SKFPTRFGGQ IRGFSSEGYI DGKNERRLDD CLKYCIVAGK KALESANLGG DKLNTIDKRK AGVLVGTGMG GLTVFSEGVQ NLIKGHRRR SPFFIPYAIT NMGSALLAID LGLMGPNYSI STACATSNYC FYAAANHIRR GEADMMIAGG TEAAIPIGL GGFVACRALS QRNDPQTAS RPWDKARDGF VMGEGAGVLV MESLEHAMKR GAPIVAEYLG GAVNCDAAHHM TDPRADGLGV SSCIERCLED AGVSPEEVNY INAHATSTLA GDLAEINAIAK KVFKSTSGIK INATKSMIGH CLGAAGGLEA IATVKAIN TG WLHPSINQFN PEQAVDFDTV PNEKKQHEVD VAISNSFGFG GHNSVVAFSA FKP
Source	Baculovirus
Target Names	KAS1
Protein Names	Recommended name: 3-oxoacyl-[acyl-carrier-protein] synthase I, chloroplastic EC= 2.3.1.41 Alternative name(s): Beta-ketoacyl-ACP synthase I Short name= KAS I
Expression Region	47-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.