



Recombinant *Borrelia burgdorferi* Flagellar assembly protein fliH (fliH)

Product Code	CSB-MP346004BUD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52611
Product Type	Recombinant Protein
Immunogen Species	<i>Borrelia burgdorferi</i> (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)
Purity	>85% (SDS-PAGE)
Sequence	MPKVLKYSSE VDNLVKFEFV EIAKPVFESL EIKEKEGKVY DIESQISNLK EELQLLRDEK LQLEELAKR QELAKEEVQI ESKRLIEEAK AKANEVLEAA KQEADLLQKE AIYKKESJET ESNAEIERLA REYEEKLKTD LEIAIAKGRE EGYSKGYESG FEDFDKVMRK LHVVIASLIA ERKGILESSS GQIVSLVMQI AIKVIKRITD SQKDIVLENV NEVLKRVKDK TQITIRVNLD DLDIVRHKKS DFISRFDIIE NLEIIEDPNI GKGGCIETN FGEIDARISS QLDKIEEKFK NFSLLS
Source	Mammalian cell
Target Names	fliH
Protein Names	Recommended name: Flagellar assembly protein fliH
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.