



Recombinant Repressor of phase 1 flagellin gene (fljA)

Product Code	CSB-YP346005SBC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52617
Product Type	Recombinant Protein
Immunogen Species	Salmonella abortus-equi
Purity	>85% (SDS-PAGE)
Sequence	MECIAVNDIS YGREAEIWPR DYSMLARRVQ FLRFNDILVR LVSNNARIIT GYIAKFNPRE NLILASDKPK GNKRIEVKLE SLAILEELSG NDAFNLSLVP ADEFNLQQYT PSRRDYFSIC NKCYKQGVGI KIYMKYGQVL TGKTTGVNAC QVGVRTSNGN HMQVMFDWVS RITSSDYAE
Source	Yeast
Target Names	fljA
Protein Names	Recommended name: Repressor of phase 1 flagellin gene
Expression Region	1-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.