



# Recombinant Solanum tuberosum Pyruvate dehydrogenase E1 component subunit alpha, mitochondrial

<b>Product Code</b>	CSB-EP346027FIG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P52903
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Solanum tuberosum (Potato)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SSDS TATITVETSL PFTSHNIDPP SRSVETSPKE LMTFFKDMTE MRRMEIAADS LYKAKLIRGF CHLYDGQEA AVGMEEAATK KDCIITAYRD HCIFLGRGGT LVEAFAELMG RRDGCSRKGG GSMHFYKKES GFYGGHGIVG AQVPLGIGLA FAQKYKKEDY VTFAMYGDGA ANQQQLFEAL NMAALWDLPA ILVCENNHYG MGTAEWRAAK SPAYYKRGDY VPGLRVDGMD VFAVKQACTF AKQHALKNGP IILEMDTYRY HGHSMSDPGS TYRTRDEISG VRQERDPVER IRSLIAHNI ATEAELKDIE KENRKVVDEA IAKAKESPMP DPSELFNTNVY VKGFGVEAYG ADRKELRATL P
<b>Source</b>	E.coli
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit alpha, mitochondrial Short name= PDHE1-A EC= 1.2.4.1
<b>Expression Region</b>	27-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.