



Recombinant *Saccharomyces cerevisiae* Cardiolipin-specific deacylase 1, mitochondrial (CLD1)

Product Code	CSB-BP346066SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53264
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	V TEKQILRTPS APTAIPLREI IYRVPSLFPR PLEDSVKDFR DFIKNEDAFQ TELLKTLPFY PTPSESKTAR LIRTVVDDEG NYINEFCIRP RKTSVPEADL KHLVFIHGYG AGLGFFIKNF EDIPLLDNEW CIHAIDLPGY GFSSRPKFPF EYPRDNIHSV QDWFHERIHT WFSKRNLLNR PEKNIVMAHS LGSYLMALYL QKYKESPSFK KLILCSPAGV SYRDFNNTAS EVEKWKPPPW WYVKLWDRNI SPFTLVRNFR QLGSKITSGW SYRRFKHILN GDPEQSKRFE ALHRYAYAIF NKRGSGEYLL SFALKCGGEP RLSLEQQLFD GKKS DILKNS NCDWLWLYGD DDWMDVNGGL RVS RFLKEKL KQKSNVIIVP HSGHHLYLDN YKFFNNILTK EMQKI
Source	Baculovirus
Target Names	CLD1
Protein Names	Recommended name: Cardiolipin-specific deacylase 1, mitochondrial EC=3.5.1.-
Expression Region	40-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.