



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YGR117C (YGR117C)

Product Code	CSB-EP346067SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53270
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MADNSTSRFY TNLIIANYLK HNGLEDTLAA FIRETALPLS ALEKSDSSNS NVGEIPLDL QSVVEDRIYY KRRSFKDRFK TLSINDDLAP IDNAKYGIQP WNHSLKFSID VKLNKSLPKD TLFISATFTE DSKYILLSSA TGYLVIYDIE KATSKSFKIN EKVKSIKLY GPIGSSGYQY VCPMNGSFYL LNPDFELVNN AVWKIHARMI THIKICNVTE SSWFVITSGM DNFLRLSLE IKNGNTFLTK LSEIKLASNC TSLNVIANGD GNGQNSFSVF LTRAEYTHIA CYSIIDAKNL VHSYNIALNN AEFSTYAFNI RDVMAVDYVH SNTKDTIGLS PSTMLVVATS HKPYMRLILV EIPMNTGHPK AMKLDKVQTY YDKILRNFAT EIYQDDFSLP ILGKLESSNG VLVGNDEGIY SVDLMTGDSR ILNIPGEANS LHDRKCMDI SKDQMRMVAG TSTKSIYILN VIRNAQ
Source	E.coli
Target Names	YGR117C
Protein Names	Recommended name: Uncharacterized protein YGR117C
Expression Region	1-476
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.