



Recombinant *Saccharomyces cerevisiae* N-acetyltransferase SLI1 (SLI1)

Product Code	CSB-EP346071SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53304
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MNLKLSAIES YFFHRSRLNL HSCFYVGIKL NELPKKSQLI AALKYTVIQH ERLTCNVFYD ELKKENFLQN ILEPLKFCDL VEYRHDWDQL GETEINHIFQ RYNFSYNENK PLWKILILPN QNQMLLLTDH VLMDGMSAIH VWETFMEGLQ MQQPVEIDET IYSPSLNSST EKIMSAPLYG DWPIPNWHI VRQLVSRRLHY WFPQTVVKNN RNLIQFANYS FPKDLLDDKP SDGTQKYKVK NTNHQWEFRL SPTHLNDILQ ECKANNTSLT SLLGALVCTS FEKIAAHEYT GSFLKIELPM NIRKPFERVL KLPSSDKLAV GNFIAVIEFN HKLHQNRGIW DIASQIQRAI RSSSEDKIID KVNEVKLLEV ISSQQYIEDK ISLNNGPSST FEVTNLGFQT FKDACNTSLP FYIVDATFNE PQGISSIFSL SVISTPGNGL HCCISYPNTL TKVLEPHWQY MKDYLNLY
Source	E.coli
Target Names	SLI1
Protein Names	Recommended name: N-acetyltransferase SLI1 EC= 2.3.1.- Alternative name(s): Sphingosine-like immunosuppressant resistant protein 1
Expression Region	1-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.