



Recombinant Danio rerio Insulin gene enhancer protein isl-2b (isl2b)

Product Code	CSB-MP346080DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53407
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MVDIIFSSSF LDDMGDHSKK KSGLAMCVGC GSQIHDQYIL RVSPDLEWHA ACLKCVECNQ YLDETCTCFV RDGKTYCKRD YVRLFGIKCA KCTLGFSSSD LVMRARDSVY HIECFRCSVC SRQLLPGDEF SVRDEELLCR ADHGLALERG PGGSPLSPGN IHTRGLHMAA DPVSVRQTPH RNHVHKQSEK TTRVRTVLNE KQLHTLRTCY NANPRPDALM KEQLVEMTGL SPRVIRVWFQ NKRCCKDKKRS ILMKQLQQQH GDKTNLQGMT GTALVAGSPI RHNPSVPGHP VDVQAYQPPW KALSEFALQS DLDQPAFQQL VSFSESGSLG NSSGSDVTSL SSQLPDT PNS MVPSPVET
Source	Mammalian cell
Target Names	isl2b
Protein Names	Recommended name: Insulin gene enhancer protein isl-2b Short name= Islet-2B Alternative name(s): Insulin gene enhancer protein isl-3 Short name= Islet-3
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.