



Recombinant *Saccharomyces cerevisiae* RNA polymerase I-specific transcription initiation factor RRN9 (RRN9)

Product Code	CSB-EP346083SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53437
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSDLDEESQI ETQIDAPIED IIRGSELTTT TADKETLKSA NELLDSLEHS HRVDLSLHLY SAYLLKRLLY KANEKKHFYE VNQFVKTKQIK DNWTSWPNPN TIIDPSVDKL YEDIPEGIAN VSVQPGEISN RALMHASDMM RVELDAQWQK FLSKSALDHD VTLDVDELNI PNEISRNILV KLDSLFEGLH DKIAKENEFD VRQDKHSNNI RANQIDDEPM QANRRIKYTY HDLVSRGCEM NEDMTDIYMK SLELYNDIPE KYKKRKFRPL KQILKKYHQP KKTSSYLKEL LSKTREDFIP VEKLLKDKRL TSKDKSKLQR LNREETEDAL NKRTFFQVKG YLEDENEISD YELDDCLIEL PNGNI
Source	E.coli
Target Names	RRN9
Protein Names	Recommended name: RNA polymerase I-specific transcription initiation factor RRN9
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.