



# Recombinant Bacillus subtilis UPF0157 protein YqkA (yqkA)

<b>Product Code</b>	CSB-BP346184BRJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P54564
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bacillus subtilis (strain 168)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIEFLQHKEA TIAREILAVQ IPAYIKETEI IGFDGIPALQ EKVKDIQASH ERFAGYRKNG KLIGVISYEK NMHYLTICRL VVHPDGFRKG IGRALLQFVI DHNEDAIEKIE VVTAENNTPA VSPLYTQAGFQ KAETVCAAQ LLLSVFHLYP KRKVEIILYN EKWAECFDEE KERLKLVFPG EIIAIHHIGS TSIPNMAAKP IIDMLIEVRS IEAVSQFDEQ MKANGYTPKG ENGIAGRRYF QKGGNKRTHH VHMYEQGNPA IERHLLFRDY LRAHPNIAKE YAVLKKQLAA QHPDSINQYI QGKDEWIKTA EENAKRWKKG RNNANGSIVC YNSENDENGG FTL
<b>Source</b>	Baculovirus
<b>Target Names</b>	yqkA
<b>Protein Names</b>	Recommended name: UPF0157 protein YqkA
<b>Expression Region</b>	1-343
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.