



Recombinant *Deinococcus radiodurans* Transcription termination factor Rho (rho)

Product Code	CSB-EP346352DJA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52153
Product Type	Recombinant Protein
Immunogen Species	<i>Deinococcus radiodurans</i> (strain ATCC 13939 / DSM 20539 / JCM 16871 / LMG 4051 / NBRC 15346 / NCIMB 9279 / R1 / VKM B-1422)
Purity	≥85% (SDS-PAGE)
Sequence	MTVTEVAPQA LPFQELQEKI LPELHLLAAG LGIENYRKLK KDALALAIME KQADAEGQSL ARGYLDITSD GYGFLQADLL DPASRSVLVT AGVIKQYHLR TGDEVIGRAR KPRENERYGS LVRVEAVNGL DPEAARQRPR FDDLTPTFPD QQLVLEDPST DDGLSLRVVD LLVPIGRGQR ALIVAPPKAG KTTLLKIAN SITKNYPDVT VMVLLVDERP EEVTDRESV QGAQVIASF DEPPQHHRV AEFVHERARR IVEEGGHVVI LLDSITRLAR ANNLVTPPTG RTLSGGLDSN ALHWPKRFLG AARNIREGGS LTI LATALVE TGSRMDDVIF EEFKGTGNAE LVLSRRLEER RIFPALDILK SGTRREELL QPEVLKMMWL LRKVISDMDP ADAMEMLLGR MGKTRNNVEF LAALAG
Source	E.coli
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.