



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YGL081W (YGL081W)

Product Code	CSB-MP346433SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53156
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGDIRTFVFA IEDTETTQGL CKTIGRSSF DQNSLCKPYN LYFDEPELSR QHAVLCIKTP IPKIEGVPSI EQLRICIRDL NNKTGTVNLV SDGPNDEIDL KNGDAFGLIA IDNHPFRDNH HLAAKLIFRI ELEYFDEARE IVKCTITNVT FGKNNTVSSF PIHSATSTED SDSSWYGLSE ASTQTEVADE CHETNTILTR GGRFSILSLR KRGSKQDQKI CSNFDRKIHE TSSFEEIEV CTDTDTTEEK EEEEKEEGD DEEGEIELEI IRVKRIKGR T KIKKTLTCFS KNKKIIPQH SNSMWLLLIV ILIFDRLLSN
Source	Mammalian cell
Target Names	YGL081W
Protein Names	Recommended name: Uncharacterized protein YGL081W
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.