



Recombinant *Saccharomyces cerevisiae* Cold sensitive U2 snRNA suppressor 2 (CUS2)

Product Code	CSB-MP346486SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53830
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDADELELKG HLKCLKKEEL LRRKQLKESN LQKRELEYNN ASKNTSIYIS GLPTDKTTKE GLTEQFCKYG MIRTNRDGEPLCKLYVNDKG AFKGDALITY SKEESVTLAI EMMNESIFLG KQIRVERAQF QNKEGDNMHG KENDLKEFNG PEPPIKRLKK AKSEGEGEVI DYNDDESLAK ADRTVIFANV FNIYKSYTND DINDIQEDLL EGCEEIGQVD SISVSPNKGE ATVVFKNNKV ALQCCKIMTG RYFDGQKLLA FISGDENTSS TSDKNEDSEV EDDL
Source	Mammalian cell
Target Names	CUS2
Protein Names	Recommended name: Cold sensitive U2 snRNA suppressor 2
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.