



# Recombinant Mouse Interferon-induced, double-stranded RNA-activated protein kinase (Eif2ak2)

<b>Product Code</b>	CSB-BP007511MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q03963
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ASDTPGFYM DKLNKYRQMH GVAITYKELS TSGPPHRRF TFQVLIDEKE FPEAKGRSKQ EARNAAKLA VDILDNENKV DCHTSASEQG LFGVNYIGLV NSFAQKKKLS VNYEQCEPNS ELPQRFICKC KIGQTMYG TG SGVTKQEAKQ LAAKEAYQKL LKSPPKTAGT SSSVVTSTFS GFSSSSSMTS NGVSQSAPGS FSSENVFTNG LGENKRKSGV KVSPDDVQRN KYTLDFARNS DFEDIEEIGL GGFGQVFKAK HRIDGKRYAI KRVKYNTEKA EHEVQALAE L NHVNIVQYHS CWEGVDYDPE HSMSDTSRYK TRCLFIQMEF CDKGTLEQWM RNRNQSKV DK ALILDLYEQI VTGVEYIHSK GLIHRDLKPG NIFLVDERHI KIGDFGLATA LENDGKSRTR RTGTLQYMSP EQLFLKHYGK EVDIFALGLI LAELLHTCFT ESEKIKFFES LRKGDFSNDI FDNKEKSLK KLLSEKPKDR PETSEILKTL AEWRNISEKK KRNTC
<b>Source</b>	Baculovirus
<b>Target Names</b>	Eif2ak2
<b>Protein Names</b>	Recommended name: Interferon-induced, double-stranded RNA-activated protein kinase EC= 2.7.11.1 Alternative name(s): Eukaryotic translation initiation factor 2-alpha kinase 2 Short name= eIF-2A protein kinase 2 Interferon-inducib
<b>Expression Region</b>	2-515
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.