



Recombinant *Xenopus laevis* Homeobox protein DLL-3 (dll3)

Product Code	CSB-BP346574XBE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54655
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MTGVYERRVA SIRASEFQSM MHHPSQDSPT LPESTATDSG YYSPGAAAGH PHHGYCSPTS ATYGKALNAY QYQYHGMNGA AANYSGKAYS DYGYGSPYHP HHQYSGAYNR VQPPSTQPEK EVSEPEVRMV NGKPKKIRKP RTIYSSFQLA ALQRRFQKTQ YLALPERAEL AASLGLTQTQ VKIWFQNKRS KIKKIMKNGE LPPEHSPSSS DPMACNSPQS PVVWEPQGSS RSLSHHPHVH SHPQPSSSSP ASSYLENSGV WYPTGTHLQN LQNHASLQHP LALASGTLY
Source	Baculovirus
Target Names	dll3
Protein Names	Recommended name: Homeobox protein DLL-3 Short name= XDLL-3
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.