



Recombinant Escherichia coli Uncharacterized protein yfjH (yfjH)

Product Code	CSB-BP346730ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52123
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MGIFKAKNPC TKNTIFTTSN TLIYGGFMIS LNDFYEQICR KRRDLAYHMS ECEWAVDTDV LEEDHPEIRI ELGRMREQFW SSEKIGTRVR LYSCDVPWET RHHTVNGQLE IKEEYTELYD PAQECWKNLS SNLTKETFLP LVIEPFSIND IFKAHLMFAS ISFFWGKSIM SENENVAFKA FHRAAELFDK CIGMTWFNIS VCNQKKLSEV RRSAGKKGGK SKAEVYHIIQ LKLVELINDS VPNDGWKNKV VAVNELIEPL WDFIQMSEFE INNQNKKYRV ATMSQDALVD TILNQWSLKN EDVKQAFDSA VRRKKRSK
Source	Baculovirus
Target Names	yfjH
Protein Names	Recommended name: Uncharacterized protein yfjH
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.