



# Recombinant Human Rho-related GTP-binding protein RhoD (RHOD)

<b>Product Code</b>	CSB-YP019688HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O00212
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTAAQAAGEE APPGVRSVKV VLVGDGGCGK TSLLMVFADG AFPESYTPTV FERYMVNLQV KGKPVHLHIW DTAGQDDYDR LRPLFYPDAS VLLLCFDVTS PNSFDNIFNR WYPEVNHFKK KVPIIVVGCK TDLCCKDKSLV NKLRRNGLEP VTYHRGQEMA RSVGAVAYLE CSARLHDNVH AVFQEAEEVA LSSRGRNFWR RITQGFC
<b>Source</b>	Yeast
<b>Target Names</b>	RHOD
<b>Protein Names</b>	Recommended name: Rho-related GTP-binding protein RhoD Alternative name(s): Rho-related protein HP1 Short name= RhoHP1
<b>Expression Region</b>	1-207
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Ras homolog, or Rho, proteins interact with protein kinases and may serve as targets for activated GTPase. They play a critical role in muscle differentiation. This protein binds GTP and is a member of the small GTPase superfamily. It is involved in endosome dynamics and reorganization of the actin cytoskeleton, and it may coordinate membrane transport with the function of the cytoskeleton.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.