



Recombinant *Saccharomyces cerevisiae* Ubiquitin fusion degradation protein 1 (UFD1)

Product Code	CSB-MP346812SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53044
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MFSGFSSFSG GNGFVNMPQT FEEFFRCYPI AMMNDRIRKD DANFGGKIFL PPSALSKLSM LNIRYPMLFK LTANETGRVT HGGVLEFIAE EGRVYLPQWM METLGIQPGS LLQISSTDVP LGQFVKLEPQ SVDFLDISDP KAVLENVLRN FSTLTVDDVI EISYNGKTFK IKILEVKPES SSKSICVIET DLVTDFAPPV GYVEPDYKAL KAQQDKEKKN SFGKGQVLDP SVLGQGSMSR RIDYAGIANS SRNKLSKFVG QGQNISGKAP KAEPKQDIKD MKITFDGEPK KLDLPEGQLF FGFPMVLPKE DEESAAGSKS SEQNFQQQGI SLRKSNNRKT KSDHDSKSK APKSPEVIEI D
Source	Mammalian cell
Target Names	UFD1
Protein Names	Recommended name: Ubiquitin fusion degradation protein 1 Short name= UB fusion protein 1 Alternative name(s): Polymerase-interacting protein 3
Expression Region	1-361
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.