



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor CWC24 (CWC24)

Product Code	CSB-BP346888SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53769
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MFRKRLVNKS SSDEKNQKKR QKINFSEEKL VASDEEKGSS DLMSLAKSGN SRTLQLSHEN EGKLQKKGED LDKYTLTVND DSTKEDLLNF ERKELAEKAK KRRPSDDNEL VLNMSGKNKR LTKQINQPTN IRTTVLMDFQ PDVCKDYKQT GYCGYGDSCK FLHSRDDFKT GWKLNQEWNA DKEDSKAVTL DLEKIPFKCT LCKEDYKSPV VTNCGHYFCG SCFAKDMKKG TKCFICHKET HGSAKVASDL QKMLNKRKS
Source	Baculovirus
Target Names	CWC24
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC24 Alternative name(s): Complexed with CEF1 protein 24
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.