



Recombinant Human Homeobox protein EMX2 (EMX2)

Product Code	CSB-MP007657HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04743
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFQPAPKRCF TIESLVAKDS PLPASRSEDP IRPAALSYAN SSPINPFLNG FHSAAAAAAG RGVYSNPDLV FAEAVSHPPN PAVPVHPVPP PHALAAHPLP SSHSPHPLFA SQQRDPSTFY PWLIHRYRYL GHRFQGNDS PESFLLHNAL ARKPKRIRTA FSPSQLLRLE HAFEKNHYVV GAERKQLAHS LSLTETQVKV WFQNRRTKFK RQKLEEEGSD SQQKKKGTHH INRWRIATKQ ASPEEIDVTS DD
Source	Mammalian cell
Target Names	EMX2
Protein Names	Recommended name: Homeobox protein EMX2 Alternative name(s): Empty spiracles homolog 2 Empty spiracles-like protein 2
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.