



Recombinant Peptidyl-prolyl cis-trans isomerase 10 (cyn-10)

Product Code	CSB-BP347085CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52017
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSVTLHTTSG DIKIELYVDD APKACENFLA LCASDYINGC IFHRNIKDFM VQTGDPTHSG KGGESIWGGP FEDEFVSALK HDSRGCVSMA NNGPDSNRSQ FFITYAKQAH LDMKYTLFGK VIDGFDTLEE IETIKVDNKY RPLVQQKIQN VTIHANPMAA D
Source	Baculovirus
Target Names	cyn-10
Protein Names	Recommended name: Peptidyl-prolyl cis-trans isomerase 10 Short name= PPIase 10 EC= 5.2.1.8 Alternative name(s): Cyclophilin-10 Rotamase 10
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.