



Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor CWC23 (CWC23)

Product Code	CSB-MP347150SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52868
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPGHELEDVI NQRLNLYDVL ELPTPLDVHT IYDDLDPQIKR KYRTLALKYH PDKHPDNPSI IHKFHLLSTA TNILTNADV R PHYDRWLIEF LRKTNDIERN KLIQKLEESE SSTIPTTTTPH PDLLQIQRHG ELLRKLKHFN LPYGDWKHLN TQDQENASQH PYYDCSTLRI VLDNFLQSN N KSNCLSHLRN QVFITLSANE IYDIYFSERN NYSKDDSI I YTVFDTPITA QHVFRNWSSG NLIPTVKDIS PLIPLHYYS D FNLETELNDD IARLVSNEPI LLD
Source	Mammalian cell
Target Names	CWC23
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC23 Alternative name(s): Complexed with CEF1 protein 23
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.