



Recombinant *Saccharomyces cerevisiae* Regulatory protein MIG2 (MIG2)

Product Code	CSB-EP347162SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53035
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPKKQTNFPV DNENRPFRC D TCHRGFHRLE HKKRHLRTH T GEKPHHCAFP GCGKSF SRSD ELKRHMRTHT GQSQRRLKKA SVQKQEF LTV SGIPTIASGV MIHQPIPV L PANMAINVQA VNGGNIHAP NAVHPMVIPI MAQPAPIHAS AASFQPATSP MPISTYTPVP SQSFTSFQSS IGSIQSNSDV SSIFSNMNV R VNTPRSVPNS PNDGYLHQQH IPQQYQHQT A SPSVAKQQKT FAHSLASALS TLQKRTPVSA PSTTIESPSS PSDSSHTSAS SSAISLPFSN APSQLAVAKE LESVYLDNR YTTKTRRERA KFEIPEEQEE DTNNSSSGSN EEEHESLDHE SSKSRKKLSG VKLPPVRNLL KQIDVFNGPK RV
Source	E.coli
Target Names	MIG2
Protein Names	Recommended name: Regulatory protein MIG2
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.