



Recombinant *Saccharomyces cerevisiae* Anaphase-promoting complex subunit DOC1 (DOC1)

Product Code	CSB-YP347165SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53068
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDPIGINKVL DHLAPSELIK PVKSCHNKPS VLVLDDRIVD AATKDLYVNG FQEEIQYQNP TPENLQHMFH QGIEILDSAR MINVTHLALW KPSSFKLGNP VDFALDDNYD TFWQSDGGQP HQLDIMFSKR MDICVMAIFF SMIADESYAP SLVKVYAGHS PSDARFYKML EVRNVNGWVA LRFLDNREDD QLLKCQFIRL LFPVNHENGK DTHLRGIRLY VPSNEPHQDT HEWAQTLPET NNVFQDAILR
Source	Yeast
Target Names	DOC1
Protein Names	Recommended name: Anaphase-promoting complex subunit DOC1 Alternative name(s): Destruction of cyclin B protein 1
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.