



# Recombinant Mouse Arginase-2, mitochondrial (Arg2)

<b>Product Code</b>	CSB-YP002006MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O08691
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VHLSVAIVG APFSRGQKKL GVEYGPAAIR EAGLLKRLSR LGCHLKDFGD LSFTNVPQDD PYNLNVVYPR SVGLANQELA EVVSRVSSGG YSCVTMGGDH SLAIGTIIGH ARHRPDLCVI WVDAHADINT PLTTVSGNIH GQPLSFLIKE LQDKVPQLPG FSWIKPCLSP PNIVYIGLRD VEPPEHFILK NYDIQYFSMR EIDRLGIQKV MEQTFDRDIG KRQRPIHLSF DIDAFDPKLA PATGTPVVG LTYREGVYIT EEIHNTGLLS ALDLVEVNPV LATSEEEAKA TARLAVDVIA SSFGQTREGG HIVYDHLPTP SSPHESENEE CVRI
<b>Source</b>	Yeast
<b>Target Names</b>	Arg2
<b>Protein Names</b>	Recommended name: Arginase-2, mitochondrial EC= 3.5.3.1 Alternative name(s): Kidney-type arginase Non-hepatic arginase Type II arginase
<b>Expression Region</b>	23-354
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Arginase catalyzes the hydrolysis of arginine to ornithine and urea. At least two isoforms of mammalian arginase exists (types I and II) which differ in their tissue distribution, subcellular localization, immunologic crossreactivity and physiologic function. The type II isoform encoded by this gene, is located in the mitochondria and expressed in extra-hepatic tissues, especially kidney. The physiologic role of this isoform is poorly understood; it is thought to play a role in nitric oxide and polyamine metabolism. Transcript variants of the type II gene resulting from the use of alternative polyadenylation sites have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.