



Recombinant *Xenopus laevis* Homeobox protein DLL-4 (dll4)

Product Code	CSB-YP347231XBE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53775
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MTGVFDSLAA DMHSNHMPSS SYHSLHKSQE SPTLPVSTAT DSSYYTNQQQ QQHCGAVSPY GQLGSYQFHG AAINGISYST KSYDLSYTG S YSSYGPYGT S PSPPHNDQEK EDCEPEVRMV NGKPKKVRKP RTIYSSFQLA ALQRRFQKTQ YLALPERAEL AASLGVTQTQ VKIWFQNRRS KFKKMWKS GE IPSDQLPVGS ESSPCSSPSA STATWDFGPH QRLQGATNGS ALQSSNSASS PSPFLSNYSW YQTSNSAPHL QSNPLLQQHH LHHHPPAPIS AGTIF
Source	Yeast
Target Names	dll4
Protein Names	Recommended name: Homeobox protein DLL-4 Short name= XDLL-4
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.