



Recombinant UPF0001 protein in piIT-proC intergenic region

Product Code	CSB-YP347292VEU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P52055
Product Type	Recombinant Protein
Immunogen Species	Vibrio alginolyticus
Purity	>85% (SDS-PAGE)
Sequence	MSSIQQNIEH ITSQIRYDEQ KCGRTPESVQ LLAVSKTKPV EAILEAYQAG QTAFGENYVQ EGVSKVQHFA EHYPDNRIEW HFIGPIQSNK SRLVAEHFDW VHTIDRTKIA QRLNDQRPSE LKPLQVLIQV NTSGEASKSG VTEAEVFELA ELISRLPNLT LRGLMSIPAN VSDYESQLHE FQKLATLKQT LEAQFPEIDT LSMGMSGDMT AAIEAGSTMV RIGTAIFGAR DYS
Source	Yeast
Protein Names	Recommended name: UPF0001 protein in piIT-proC intergenic region Alternative name(s): ORF1
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.