



Recombinant *Saccharomyces cerevisiae* Protein PXR1 (PXR1)

Product Code	CSB-EP347407SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53335
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGLAATRTKQ RFGLDPRNTA WSNDTSRFGH QFLEKFGWKP GMGLGLSPMN SNTSHIKVSI KDDNVGLGAK LKRKDKKDEF DNCECAGLDV FQRILGRLNG KESKISEELD TQRKQKIIDG KWGIHFVKGE VLASTWDPKT HKLRNYSNAK KRKREGDDSE DEDDDDKEDK DSDKKKHKKH KKHKKDKKKD KKDKKEHKKH KKEEKRLKKE KRAEKTETK KTSKLSSES ASNIPDAVNT RLSVRSKWIK QKRAALMDSK ALNEIFMITN D
Source	E.coli
Target Names	PXR1
Protein Names	Recommended name: Protein PXR1 Alternative name(s): G-patch nucleolar protein PinX1-related protein 1
Expression Region	1-271
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.