



Recombinant *Saccharomyces cerevisiae* Multifunctional methyltransferase subunit TRM112 (TRM112)

Product Code	CSB-EP347430SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53738
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKFLTTNFLK CSVKACDTSN DNFPLQYDGS KCQLVQDESI EFNPEFLLNI VDRVDWPAVL TVAAELGNNALPPTKPSFPS SIQELTDDDM AILNDLHTLL LQTSIAEGEM KCRNCGHIYY IKNGIPNLLL PPHLV
Source	<i>E.coli</i>
Target Names	TRM112
Protein Names	Recommended name: Multifunctional methyltransferase subunit TRM112 Alternative name(s): eRF1 methyltransferase subunit TRM112 Short name= eRF1 MTase subunit TRM112 tRNA methyltransferase 112
Expression Region	1-135
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.