



Recombinant *Saccharomyces cerevisiae* Adenine deaminase (AAH1)

Product Code	CSB-EP347445SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53909
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVSVEFLQEL PKCEHHLHLE GTLEPDLLFP LAKRNDIILP EGFPKSVEEL NEKYKKFRDL QDFLDYYYIG TNVLISEQDF FDLAWAYFKK VHKQGLVHAE VFYDPQSHTS RGISIETVTK GFQRACDKAF SEFGITSKLI MCLLRHIEPE ECLKTIEEAT PFIKDG TISA LGLDSA EKPFP PHLFVECYG KAASLNKDLK LTAHAGEEGP AQFVSDALDL LQVTRIDHGI NSQYDEELLD RLSRDQTMLT ICPLSNVKLQ VVQSVSELPL QKFLDRDVPF SLNSDDPAYF GGYILDVYTQ VSKDFPHWDH ETWGRIAKNA IKGSWCDDKR KNGLLSRVDE VVTKYSH
Source	E.coli
Target Names	AAH1
Protein Names	Recommended name: Adenine deaminase Short name= ADE EC= 3.5.4.2 Alternative name(s): Adenine aminohydrolase Short name= AAH
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.