



# Recombinant Human Eukaryotic translation initiation factor 1A, Y-chromosomal (EIF1AY)

<b>Product Code</b>	CSB-YP007507HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O14602
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	PKNKGKGGK NRRRGKNENE SEKRELVFKE DGQEYAQVIK MLGNRLEAL CFDGVKRLCH IRGKLRKKVW INTSDIILVG LRDYQDNKAD VILKYNAD RSLKAYGELP EHAKINETDT FGPGDDDEIQ FDDIGDDDED IDDI
<b>Source</b>	Yeast
<b>Target Names</b>	EIF1AY
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 1A, Y-chromosomal Short name= eIF-1A Y isoform Alternative name(s): Eukaryotic translation initiation factor 4C Short name= eIF-4C
<b>Expression Region</b>	2-144
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a protein similar to eukaryotic translation initiation factor 1A (EIF1A). EIF1A is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNA <sub>i</sub> and eIF3) to the 5' end of capped RNA.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.