



Recombinant Human Actin-binding LIM protein 1 (ABLIM1)

Product Code	CSB-BP001107HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14639
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MPAFLGLKCL GKLCSSSEKSK VTSSERTSAR GSNRKRLIVE DRRVSGTSFT AHRRTATHL LYLCPKDYCP RGRVCNSVDP FVAHPQDPHH PSEKPIVHCH KCGEPCKGEV LRVQTKHFHI KCFTCKVCGC DLAQGGFFIK NGEYLCTLDY QRMYGTRCHG CGEFVEGEVV TALGKTYHPN CFACTICKRP FPPGDRVTFN GRDCLCQLCA QPMSSSPKET TFSSNCAGCG RDIKNGQALL ALDKQWHLGC FKCKSCGKVL TGEYISKDGA PYCEKDYQGL FGVKCEACHQ FITGKVLEAG DKHYHPSCAR CSRCNQMFTE GEEMYLQGST VWHPDCKQST KTEEKLRPTR TSSESIYSRP GSSIPGSPGH TIYAKVDNEI LDYKDLAIP KVKAIYDIER PDLITYEPFY TSGYDDKQER QSLGESPRTL SPTPSAEGYQ DVRDRMIHRS TSQGSINSPV YSRHSYTPPT SRSPQHFRP GNEPSSGRNS PLPYRPDSRP LTPTYAQAPK HFHVPDQGIN IYRKPIYKQ HAALAAQSKS SEDIKFSKF PAAQAPDPSE TPKIETHWP GPPSFAVVG DMKRRSSGRE EDDEELLRRR QLQEEQLMKL NSGLGQLILK EEMEKESSRER SLLASRYDS PINSASHIPS SKTASLPGYG RNGLHRPVST DFAQYNSYGD VSGGVRDYQT LPDGHMPAMR MDRGVSMPTM LEPKIFPYEM LMVTNRGRNK ILREVDTRL ERHLAPEVFR EIFGMSIQEF DRLPLWRRND MKKKAKLF</p>
Source	Baculovirus
Target Names	ABLIM1
Protein Names	<p>Recommended name: Actin-binding LIM protein 1 Short name= abLIM-1 Alternative name(s): Actin-binding LIM protein family member 1 Actin-binding double zinc finger protein LIMAB1 Limatin</p>
Expression Region	1-778
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	<p>This gene encodes a cytoskeletal LIM protein that binds to actin filaments via a domain that is homologous to erythrocyte dematin. LIM domains, found in over 60 proteins, play key roles in the regulation of developmental pathways. LIM domains also function as protein-binding interfaces, mediating specific protein-protein interactions. This protein could mediate such interactions between actin filaments and cytoplasmic targets. Alternatively spliced transcript variants</p>



encoding different isoforms have been identified.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.