



Recombinant Human Inositol monophosphatase 2 (IMPA2)

Product Code	CSB-MP011695HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14732
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKPSGEDQAA LAAGPWEECF QAAVQLALRA GQIIRKALTE EKRVTSTKTS ADLVTETDHL VEDLISEL RERFPSHRFIA EAAAASGAKC VLTHSPTWII DPIDGTCNFV HRFPTVAVSI GFAVRQELEF GVIYHCTEER LYTGRRRGRGA FCNGQRLRVS GETDLSKALV LTEIGPKRDP ATCLKFLSNM ERLHAKAHG VRVIGSSTLA LCHLASGAAD AYYQFGLHCW DLAAATVIIR EAGGIVIDTS GGPLDLMACR VVAASTREMA MLIAQALQTI NYGRDDEK
Source	Mammalian cell
Target Names	IMPA2
Protein Names	Recommended name: Inositol monophosphatase 2 Short name= IMP 2 Short name= IMPase 2 EC= 3.1.3.25 Alternative name(s): Inositol-1(or 4)-monophosphatase 2 Myo-inositol monophosphatase A2
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.