



Recombinant Avian sarcoma virus Transforming protein Qin (V-QIN)

Product Code	CSB-BP348043ATF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56260
Product Type	Recombinant Protein
Immunogen Species	Avian sarcoma virus (strain 31) (ASV31)
Purity	>85% (SDS-PAGE)
Sequence	MLDMGDRKEV KMLPKSSFSI NNLVPEAVQS DNHSGHSHHN SHHPHHHHHH HHHHPPPPQQ PQRAAAAEEE DEEKAPLLLP PPAAGALEAA KAEALAGKGE AGAAAAELEE KEKAAEEKKG AAEGGKDGES GKEGEKKNKG YEKPPFSYNA LIMMAIRQSP EKRLTLNGIY EFIMKNFPYY RENKQGWQNS IRHNLSLNKC FVKVPRHYDD PGKGNWMLD PSSDDVFIGG TTDKLRRRST TSRAKLAFKR GARLTSTGLT FMDRAGSLYW PMSPFSLHH PRASSTLSYN GTASAYPSHP MPYSSVLTQN SLGNNHSFST SNGLSVDRLV NGEIPYATHH LTAALAASV PCGLSVPCSG TYSLNPCSVN LLAGQTSYFS PTSLTPQ
Source	Baculovirus
Target Names	V-QIN
Protein Names	Recommended name: Transforming protein Qin Alternative name(s): Oncogene Qin
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.