



Recombinant Human High mobility group protein B3 (HMGB3)

Product Code	CSB-MP010561HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15347
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAKGDPKKPK GKMSAYAFFV QTCREEHKKK NPEVPVNFAE FSKKCSERWK TMSGKEKSKF DEMAKADKVR YDREMKDYGP AKGGKKKKDP NAPKRPPSGF FLFCSEFRPK IKSTNPGISI GDVAKKLGEM WNNLNDSEKQ PYITKAAKLK EKYEKDVADY KSKGKFDGAK GPAKVARKKV EEEEEEEEEE EEEEEEEEDE
Source	Mammalian cell
Target Names	HMGB3
Protein Names	Recommended name: High mobility group protein B3 Alternative name(s): High mobility group protein 2a Short name= HMG-2a High mobility group protein 4 Short name= HMG-4
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.