



# Recombinant Human Eukaryotic translation initiation factor 3 subunit H (EIF3H)

<b>Product Code</b>	CSB-EP007537HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O15372
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASRKEGTGS TATSSSSTAG AAGKGKGGKGG SGDSAVKQVQ IDGLVVLKII KHYQEEGQGT EVVQGVLLGL VVEDRLEITN CFPFPQHTED DADFDEVQYQ MEMMRS LRHV NIDHLHVGWY QSTYYGSFVT RALLDSQFSY QHAIEESVVL IYDPIKTAQG SLSLKAYRLT PKLMEVCKEK DFSPEALKKA NITFEYMFEE VPIVIKNSHL INVLMWELEK KSAVADKHEL LSLASSNHLG KNLQLLMDRV DEMSQDIVKY NTYMRNTSKQ QQQKHQYQQR RQQENMQRQS RGEPLPEED LSKLFKPPQP PARMDSLIA GQINTYCQNI KEFTAQNLGK LFMAQALQEY NN
<b>Source</b>	E.coli
<b>Target Names</b>	EIF3H
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3 eIF-3-gamma eIF3 p40 subunit
<b>Expression Region</b>	1-352
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.