



Recombinant Human TOX high mobility group box family member 3 (TOX3)

Product Code	CSB-BP024075HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15405
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDVRFYPAAA GDPASLDFAQ CLGYYGYSKF GNNNNYMNMA EANNAFFAAS EQTFHTPSLG DEEFEIPPIT PPESDPALG MPDVLLPFQA LSDPLPSQGS EFTPQFPPQS LDLPSITISR NLVEQDGV LH SSSLHMDQSH TQVSQYRQDP SLIMRSIVHM TDAARSGVMP PAQLTTINQS QLSAQLGLNL GGASMPHTSP SPPASKSATP SPSSSINEED ADEANRAIGE KRAAPDSGKK PKTPKKKKKK DPNEPQKPVS AYALFFRDTQ AAIKGQNPNA TFGEVSKIVA SMWDSLGE EQ KQVYKRKTEA AKKEYLKALA AYRASLVSKA AAESAEAQTI RSVQQTLAST NLTSSLLLNT PLSQHGTVSA SPQTLQQSLP RSIAPKPLTM RLPMNQIVTS VTIAANMPSN IGAPLISSMG TTMVGSAPST QVSPSVQTQQ HQMQLQQQQQ QQQQQMQMQQ QQQLQQHQM H QQIQQMQMQQ HFQHMQQHL QQQQHLQQQ INQQQLQQQL QQRLQLQQQL HMQHQSQPSP RQHSPVASQI TSPIPAIGSP QPASQQHQSQ IQSQTQTQVL SQVSIF
Source	Baculovirus
Target Names	TOX3
Protein Names	Recommended name: TOX high mobility group box family member 3 Alternative name(s): CAG trinucleotide repeat-containing gene F9 protein Trinucleotide repeat-containing gene 9 protein
Expression Region	1-576
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.