



Recombinant Solanum tuberosum Aspartic protease inhibitor 9

Product Code	CSB-EP348261FIG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P58521
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	ESPLPKPVLD TNGKELNPNS SYRIISIGAG ALGGDVYLGK SPNSDAPCPD GVFRYNSDVG PSGTPVRFIP LSGGIFEDQL LNIQFNIPTV KLCVSYTIWK VGNLNAYFRT MLLETGGTIG QADNSYFKIV KLSNFGYNLL SCPFTSIICL RCPEDQFCAK VGVVIQNGKR RLALVNNENPL DVLFQEV
Source	E.coli
Protein Names	Recommended name: Aspartic protease inhibitor 9 Alternative name(s): Novel inhibitor of cathepsin D Short name= NID
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.