



Recombinant Acetylornithine deacetylase (argE)

Product Code	CSB-MP348666SZB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P59600
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MKNKLPPFIE IYRALIATPS ISATEEALDQ SNADLITLLA DWFKDLGFNV EVQPVPGTRN KFNMLASYGQ GAGGLLLAGH TDTVPFDDGR WTRDPFTLTE HDGKLYGLGT ADMKGGFAFI LDALRDVDVT KLKPLYILA TADEETSMAG ARYFAETTAL RPDCAIIGEP TSLQPVRAHK GHISNAIRIQ GQSGHSSDPA RGINAIELMH DAIGHILQLR DNLKERYHYE AFTVPYPTLN LGHIHGGDAS NRICACCELH MDIRPQPGMT LNELNGLLND ALAPVSERWP GRLTVDELHP PIPGYECPPN HQLVEVVEKL LGAKTEVVNY CTEAPFIQTL CPTLVLGPGS INQAHQPDEY LETRFIKPTR ELIIQVIHHF CWH
Source	Mammalian cell
Target Names	argE
Protein Names	Recommended name: Acetylornithine deacetylase Short name= AO Short name= Acetylornithinase EC= 3.5.1.16 Alternative name(s): N-acetylornithinase Short name= NAO
Expression Region	1-383
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.