



Recombinant Phalaris aquatica Major pollen allergen Pha a 5.1

Product Code	CSB-BP348695PGU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56164
Product Type	Recombinant Protein
Immunogen Species	Phalaris aquatica (Canary grass)
Purity	>85% (SDS-PAGE)
Sequence	TPPTPRT PPLLPPRAR DKATLTSRSV EDINAASRRP WWASVPPADK FKTFADHVLV VPNAVDTSA TKAPQLKAKL DAAYRVAYEA AEGSTPEAKY DAFIAALTEA LRVIAGAFEV HAVKPATEEV VADPVGELQI VDKIDAAFKI AATAANSAPA NDKFTVFEGA FNKAIKESTA GAYETYKFIP SLEAAVKQAY GATVARAPEV KYAVFEAGLT KAITAMSEAQ KVAKPPLSPQ PPQVLPLAAG GAATVAAASD VRVCRSHGTL QDACLLRCRG GCQPVVWRGG SHRARGGYKV
Source	Baculovirus
Protein Names	Recommended name: Major pollen allergen Pha a 5.1 Alternative name(s): Allergen Pha a 5 Allergen= Pha a 5.1
Expression Region	24-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.