



Recombinant Magnaporthe oryzae Scytalone dehydratase (SDH1)

Product Code	CSB-EP348697MPM-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56221
Product Type	Recombinant Protein
Immunogen Species	Magnaporthe oryzae (strain 70-15 / ATCC MYA-4617 / FGSC 8958) (Rice blast fungus) (Pyricularia oryzae)
Purity	>85% (SDS-PAGE)
Sequence	MGSQVQKSDE ITFSDYLGLM TCVYEWADSY DSKDWDRLRK VIAPTLRIDY RSFLDKLWEA MPAEEFVGMV SSKQVLGDPT LRTQHFIGGT RWEKVSEDEV IGYHQLRVPH QRYKDTTMKE VTMKGAHSA NLHWYKKIDG VWKFAGLKP D IRWGEFDFDR IFEDGRETFG DK
Source	E.coli
Target Names	SDH1
Protein Names	Recommended name: Scytalone dehydratase EC= 4.2.1.94
Expression Region	1-172
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.