



# Recombinant Mouse Bcl-2-related ovarian killer protein (Bok)

<b>Product Code</b>	CSB-YP002773MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O35425
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVLRSSVF AAEIMDAFDR SPTDKELVAQ AKALGREYVH ARLLRAGLSW SAPERASPAP GGRLAEVCTV LLRLGDELEQ IRPSVYRNVA RQLHIPLQSE PVVTDAFLAV AGHIFSAGIT WGKVVSLYSV AAGLAVDCVR QAQPAMVHAL VDCLGEFVRK TLATWLRRRG GWTDVLKCVV STDPGFRSHW LVATLCSFGR FLKAAFFLLL PER
<b>Source</b>	Yeast
<b>Target Names</b>	Bok
<b>Protein Names</b>	Recommended name: Bcl-2-related ovarian killer protein Alternative name(s): Apoptosis activator Mtd Protein matador
<b>Expression Region</b>	1-213
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This protein belongs to the BCL-2 protein family. BCL-2 family members form hetero- or homodimers and act as anti- and pro-apoptotic regulators that are involved in a wide variety of cellular activities. This protein contains all four BCL-2 like domains (BH1, 2, 3 and 4) and is a pro-apoptotic BCL-2 protein identified in the ovary.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.