



Recombinant Human Protein kinase C delta type (PRKCD)

Product Code	CSB-BP018702HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05655
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MAPFLRIAFN SYELGSLQAE DEANQPFCAV KMKEALSTER GKTLVQKKPT MYPEWKSTFD AHIEYGRVIQ IVLMRAAEEP VSEVTVGVSV LAERCKKNG KAEFWLDLQP QAKVLMSVQY FLEDVDCKQS MRSEDEAKFP TMNRRGAIKQ AKIHYIKNHE FIATFFGQPT FCSVCKDFVW GLNKQGYKCR QCNAAIHKKC IDKIIGRCTG TAANSRDTIF QKERFNIDMP HRFKVHNYMS PTFCDHCGSL LWGLVKQGLK CEDCGMNVHH KCREKVANLC GINQKLLAEA LNQVTQRASR RSDSASSEPV GIYQGFEEKT GVAGEDMQDN SGTYGKIWEG SSKCNINNF FHKVLGKGSF GKVLLGELKG RGEYFAIKAL KKDVLIDDD VECTMVEKRV LTLAAENPFL THLICTFQTK DHLFFVMEFL NGGDLMYHIQ DKGRFELYRA TFYAAEIMCG LQFLHSGKII YRDLKLDNVL LDRDGHKIA DFGMCKENIF GESRASTFCG TPDYIAPEIL QGLKYTFSVD WWSFGVLLYE MLIGQSPFHG DDEDELFESE RVDTPHYPRW ITKESKDILE KLFEREPTKR LGVTGNIKIH PFFKTINWTL LEKRRLEPPF RPKVKSPRDY SNFDQEFLNE KARLSYSDKN LIDSMDQSAF AGFSFVNPKE EHLLED</p>
Source	Baculovirus
Target Names	PRKCD
Protein Names	Recommended name: Protein kinase C delta type EC= 2.7.11.13 Alternative name(s): Tyrosine-protein kinase PRKCD EC= 2.7.10.2 nPKC-delta
Expression Region	1-676
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	<p>Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play distinct roles in cells. This protein is one of the PKC family members. Studies both in human and mice demonstrate that this kinase is involved in B cell signaling and in the regulation of growth, apoptosis,</p>



and differentiation of a variety of cell types. Alternatively spliced transcript variants encoding the same protein have been observed.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.