



Recombinant Arabidopsis thaliana Acetyl-coenzyme A carboxylase carboxyl transferase subunit beta, chloroplastic (accD)

| | |
|--------------------------|--|
| Product Code | CSB-MP348922DOA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P56765 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MEKSWFNFMF SKGELEYRGE LSKAMDSFAP GEKTTISQDR FIYDMDKNFY GWDERSYSS SYSNNVDLLV SSKDIRNFIS DDTFFVRDSN KNSYSIFFDK KKKIFEIDND FSDLEKFFYS YCSSSYLNNR SKGDNDLHYD PYIKDTKYNC TNHINSCIDS YFRSYICIDN NFLIDSNNFN ESYIYNFICS ESGKIRESKN YKIRTNRNRS NLISSKDFDI TQNYNQLWIQ CDNCYGLMYK KVKMNVCEQC GHYLMSSSE RIELLIDPGT WNPMDDEMVS ADPIKFHSKE EPYKNRIDSA QKTTGLTDAV QTGTGQLNGI PVALGVMDFR FMGGSMGSV GEKITRLIEY ATNQCLPLIL VCSSGGARMQ EGSLSLMQMA KISSVLCDYQ SSKKLFYISI LTSPTTGGVT ASFGMLGDII IAEPYAYIAF AGKRVIEQTL KKAVPEGSQA AESLLRKGLL DAIVPRNLLK GVLSELFQLH AFFPLNTN |
| Source | Mammalian cell |
| Target Names | accD |
| Protein Names | Recommended name: Acetyl-coenzyme A carboxylase carboxyl transferase subunit beta, chloroplastic Short name= ACCase subunit beta Short name= Acetyl-CoA carboxylase carboxyltransferase subunit beta EC= 6.4.1.2 |
| Expression Region | 1-488 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, |



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.