



Recombinant Buchnera aphidicola subsp. Acyrthosiphon pisum Transcription termination factor Rho (rho)

Product Code	CSB-EP349010BWZ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P57652
Product Type	Recombinant Protein
Immunogen Species	Buchnera aphidicola subsp. Acyrthosiphon pisum (strain APS) (Acyrthosiphon pisum symbiotic bacterium)
Purity	>85% (SDS-PAGE)
Sequence	MNLTALKNMP VSELITLGEK MGLENLARMR KQDIIFAILK QHAKSGEDIF GDGVLEILQD GFGFLRSADS SYLAGPDDIY VSPSQIRRFN LRTGDTIAGK IRPPKEGERY FALLKVNEVN YDKPENARSK ILFENLTPLH ANSRLRMERG NGSTEDLTAR VLDLAPIGR GQRGLIVAPP KAGKTILLQN IAQSIAYNHP DCVLMVLLID ERPEEVTEMQ RLVKGEVVAS TFDEPASRHV QVAEMVIEKA KRLVEHKKDV IILLDSITRL ARAYNTVPA SGKVLTGGVD ANALHRPKRF FGAARNVEEG GSLTIIATAL VDTGSKMDEV IYEEFKGTGN MELPLSRKIA EKRVFPAIDY NRSCTRKEEL LTLPDELQKM WILRKIIHPM SEIDAMEFLL NKLAMTKTND EFFDMMKRS
Source	E.coli
Target Names	rho
Protein Names	Recommended name: Transcription termination factor Rho EC= 3.6.4.- Alternative name(s): ATP-dependent helicase Rho
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.