



Recombinant UPF0012 hydrolase ybeM (ybeM)

Product Code	CSB-MP349042EOD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P58054
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	>85% (SDS-PAGE)
Sequence	MLVAAGQFAV TSVWEKNAEI CASLMAQAAE NDVSLFVLPE ALLARDDHDA DLSVKSAQLL EGEFLGRLRR ESKRNMMTTI LTIHVPSTPG RAWNMLVALQ AGNIVARYAK LHLYDAFAIQ ESRRVDAGNE IAPLLEVEGM KVGLMTCYDL RFPELALAQA LQGAEILVLP AAWVRGPLKE HHWSTLLAAR ALDTTCYMVA AGECGNKNIG QSRIIDPFGV TIAAASEMPA LIMAEVTPER VRQVRAQLPV LNNRRFAPPQ LL
Source	Mammalian cell
Target Names	ybeM
Protein Names	Recommended name: UPF0012 hydrolase ybeM EC= 3.5.-.-
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.